









Company	CareFusion		Circadiance	DeVilbiss Healthcare	
Product Name	 Advantage Series Nasal Mask	 SNAPP 2.0 Direct Nasal Mask	 SleepWeaver	 EasyFit Full Face Mask	 EasyFit Nasal Mask
Website	www.carefusion.com	www.carefusion.com	www.circadiance.com	www.devilbisshealthcare.com	www.devilbisshealthcare.com
Price/Warranty	Masks and headgear free from defects in materials and workmanship for 90 days.	Masks and headgear free from defects in materials and workmanship for 90 days.	\$170 retail/30-day money-back guarantee	30-day warranty	30-day warranty
Usage	CPAP	Yes	Yes	Yes	Yes
	APAP	Yes	Yes	Yes	Yes
	Bi-level	Yes	Yes	Yes	Yes
Size Options	S, M, L	S, M, L, XL	Single size	S, M, L, XL	S, M, L, XL
Latex Free	Yes	Yes	Yes	Yes	Yes
Intentional Leak (ml/s)	30 l/min @ 10 cm H <sub>2</sub> O	30 l/min @ 10 cm H <sub>2</sub> O	25 l/min @ 10 cm H <sub>2</sub> O	20 lpm @ 8 cm H <sub>2</sub> O, 40 lpm @ 20 cm H <sub>2</sub> O, 50 lpm @ 25 cm H <sub>2</sub> O	20 lpm @ 8 cm H <sub>2</sub> O, 40 lpm @ 25 cm H <sub>2</sub> O
Additional Features	Soft triple cushion for user comfort. Bellow notch at bridge of mask to allow for users who overtighten headgear. Stable forehead cushion distributes pressure.	Bellows nasal interface cushions for patient comfort. Mask cushion flanges provide mask stability and maintain nasal seal for restless sleepers. Conveniently packaged in kit that includes 4 sized cushions.	SleepWeaver lets you fit your patients quickly and comfortably to attain a reliable reading during titration.	Replaceable cushions and forehead pads available in gel or silicone. Exhalation vent is ultra-quiet and designed to direct air parallel to the patient's face and not at a bed partner.	Replaceable cushions and forehead pads available in gel or silicone. Exhalation vent is quiet and designed to direct air parallel to the patient's face and not at a bed partner.
Fitting Support Solutions	Sizing gauges available by request. Fitting instructions on mask insert instructions for use.	Fitting gauge included in package along with illustrations for fitting tips.	Call and order DVD or view instructional video at www.sleepweaver.com.	Headgear accommodates an array of head profiles. Adjustable forehead support with wider than average forehead cushions for improved customization.	Headgear accommodates an array of head profiles. Adjustable forehead support with wider than average forehead cushions for improved customization.
Pressure Range	4-20 cm H <sub>2</sub> O	4-20 cm H <sub>2</sub> O	1-20 cm H <sub>2</sub> O	4-30 cm H <sub>2</sub> O	4-30 cm H <sub>2</sub> O
Cleaning Instructions	Wipe cushions daily with clean wet cloth. Wash weekly, immersing mask and headgear into mild dishwashing soap solution. Dry before use.	Wipe cushions daily with clean wet cloth. Wash weekly, immersing mask and headgear into mild dishwashing soap solution. Dry before use.	Use cold or warm water. Can be washed in clothes washer with detergent or by hand with laundry detergent (always use oil-free products). Hang to dry.	Disassemble mask before cleaning. Immerse all components in warm water. Use a mild detergent and soft cloth to remove residue. Clean cavities with a soft brush. Rinse with clean water and allow to air dry.	Disassemble mask before cleaning. Immerse all components in warm water. Use a mild detergent and soft cloth to remove residue. Clean cavities with a soft brush. Rinse with clean water and allow to air dry.
Headgear Description	4-point adjustable headgear. Made from AccuMed Breath-O-Prene fabric with SilverTec lining to prohibit bacterial growth.	2-point adjustable headgear. Made from AccuMed Breath-O-Prene fabric with SilverTec lining to prohibit bacterial growth.	Headgear is made of flexible Velcro and elastic straps that allow the SleepWeaver mask to be quickly sized to any face.	The four-point headgear crosses at the back of the patient's head. In contrast to conventional headgear, HEADstrap ensures that no pressure is applied to the patient's neck. It is suitable for every head size.	The four-point headgear crosses at the back of the patient's head. In contrast to conventional headgear, HEADstrap ensures that no pressure is applied to the patient's neck. It is suitable for every head size and is easy to use.
Interface Description	Advantage Series Nasal Masks are lightweight and have a low profile. Exhalation port directs and diffuses downward, away from bed partner. Large headgear clips for uncomplicated application and disconnect. Lightweight mask with low patient profile to enhance user comfort.	SNAPP 2.0 can be worn with tubing down or up overhead. Nasal cushions form comfortable seal at nares. Quiet exhalation port directs airflow down and away from bed partner. Headgear has 3 adjustment selections, customizing the fit to individual user preference.	The SleepWeaver achieves a good seal and quiet operation using the properties of a balloon. When a balloon is inflated, all points inside the balloon are at the same pressure. This allows the surface areas of the mask to conform to the face, like a balloon, creating the seal.	The EasyFit Full Face Mask was designed by utilizing computer analyses of facial characteristics from 50,000 individuals globally. Identifying consistencies between most of these facial profiles contributed to the development of the scallop-shaped cushion. This, combined with the lightweight and durable nasal mask, offers a superior fit, seal, and comfort. It also includes the exclusive ball-and-socket joint at the tubing connector, allowing greater freedom of movement, even for the most active sleepers.	The EasyFit Nasal Mask was designed by utilizing computer analyses of facial characteristics from 50,000 individuals globally. Identifying consistencies between most of these facial profiles contributed to the development of the scallop-shaped cushion. The mask includes a ball-and-socket joint at the tubing connector, allowing greater freedom of movement.

Company		evo Medical Solutions	Fisher & Paykel Healthcare
Product Name			
	ComfortFit Full Face Deluxe Mask (100FD Series)	ComfortFit Nasal Deluxe+ Mask (100ND+ Series)	Opus <sup>360</sup> Nasal Pillows Mask
Website		www.evomedical.com	www.fphcare.com
Price/Warranty		N/A	N/A
Usage	CPAP	Yes	Yes
	APAP	Yes	Yes
	Bi-level	Yes	Yes
Size Options		S, M, L	S, M, L
Latex Free		Yes	Yes
Intentional Leak (ml/s)		N/A	30 l/min at 9 cm of pressure
Additional Features		One-piece molded silicone rubber; small, light, and comfortable. Clear mask and frame allow for a good view of the nasal area.	One-piece molded silicone rubber; small, light, and comfortable. Clear mask and frame allow for a good view of the nasal area.
Fitting Support Solutions		Includes step-by-step instructions that help the patient assemble the mask for maximum fit and comfort.	Enclosed illustrated fitting guide, fully assembled, fit-and-forget technology, Web site video clips.
Pressure Range		N/A	4-20 cm H <sub>2</sub> O
Cleaning Instructions		Gently hand wash mask with mild soap and warm water. Rinse mask and components thoroughly. Air dry before reusing.	Rinse thoroughly with fresh water; ensure all soap is removed. Dry in direct sunlight before reassembling. Weekly: Hand wash the headgear in pure soap dissolved in lukewarm water.
Headgear Description		Adjustable, soft headgear with fasteners. Extra buckles provide maximum versatility and movement.	Easy one-step headgear removal for cleaning enables reassembly without readjustments. Simply slide the Side Strap down and out of the Hidden Strap Release hooks.
Interface Description		Lightweight, soft, and comfortable. The 100FD Series offers replacement cushions, maximizing the life of the mask.	The Opus <sup>360</sup> Nasal Pillows Mask is small, quiet, and comfortable. The slim-line design allows for almost any sleeping position without mask displacement. Its 360-degree pivoting ball-and-socket elbow provides greater versatility. With its new silicone pillows and facially contoured mask frame for extra stability, patient comfort is paramount.

**Receive your own  
FREE  
SUBSCRIPTION  
to *Sleep Review*  
magazine.**



*Sleep Review* is the leading magazine dedicated exclusively to the sleep therapy and diagnostics market, and it provides business-management strategies and new product/technology development.






**Subscribe now  
or renew online at:  
[www.submag.com/sub/  
sn?pk=hse](http://www.submag.com/sub/sn?pk=hse)**






**SCAN ME!**



Scan the QR code above with your smartphone to connect to our Online Buyer's Guide.

**Sleep Review**  
THE JOURNAL FOR SLEEP SPECIALISTS

Company		Fisher & Paykel Healthcare	Hans Rudolph Inc	Invacare Corporation	Philips Respironics	
Product Name		 Zest Nasal Mask	 V2 Full Face CPAP Mask	 7500 Series NIV V2 Mask	 Twilight Full Face Mask	 EasyLife
Website		www.fphcare.com	www.rudolphkc.com	www.rudolphkc.com	www.invacare.com	www.philips.com/easylife
Price/Warranty		N/A	6-month warranty and size guarantee	\$215; 12-month warranty	90-day warranty	90-day warranty for defects of workmanship and materials
Usage	CPAP	Yes	Yes	Yes	Yes	Yes
	APAP	Yes	Yes	Yes	Yes	Yes
	Bi-level	Yes	Yes	Yes	Yes	Yes
Size Options		Petite, Standard, Plus	XS, Petite, S, M, L	P, XS, S, M, L	M	P, S, M, MW, L
Latex Free		Yes	Yes	Yes	Yes	Yes
Intentional Leak (ml/s)		32 l/min at 9 cm of pressure	N/A	N/A	N/A	4-30 standard liters per minute depending on PAP pressure.
Additional Features		With Foam Cushion FlexiFit Technology, automatically adjusts to a wide range of nasal and facial contours.	Detachable swivel port, airtight seal, new headgear.	Fits under chin, single piece soft silicone rubber, foam strap type headgear, 5 sizes, leak free.	Unlimited position of forehead support for customized fit, uses headgear as forehead cushion.	FitPack and DuoPack options
Fitting Support Solutions		Enclosed illustrated fitting guide, fully assembled, fit-and-forget technology, Web site video clips.	Call or e-mail HRI for pictured fitting instructions.	Mask sizing gauge available. Call or e-mail HRI for pictured fitting instructions.	Adjustable headgear, 360-degree swivel, adjustable forehead support.	Fitting guide, stand-alone sizing gauge, integrated sizing gauge on packaging, Web site video clips.
Pressure Range		4-20 cm H <sub>2</sub> O	3-40 cm H <sub>2</sub> O	0-40 cm H <sub>2</sub> O	3-20 cm H <sub>2</sub> O	4-30 cm H <sub>2</sub> O
Cleaning Instructions		Rinse thoroughly with fresh water; ensure all soap is removed. Dry in direct sunlight before reassembling. Weekly: Hand wash the headgear in pure soap dissolved in lukewarm water.	Mask, swivel dishwasher safe. Can hand-clean headgear. V2: steam autoclaved, liquid chemical, hot water pasteurized in a hospital setting.	Minimal simple two piece disassembly. Dishwasher safe; hot or cold water and soap. Steam autoclaved, hot water pasteurized, liquid chemicals.	Clean regularly with soap and water.	Clean with mild soapy water without perfumes or lotions. Multiple-patient use lab reprocessing instructions also available.
Headgear Description		Stretchgear Headgear Deluxe mask suspension system. A comfortable and breathable elastic headgear, specially designed to maintain a firm seal.	Simple Strap Foam Headgear with side adjustments and top of the head Tri-Glide adjustment.	3 sizes, soft, lightweight, simple strap type foam headgear.	Breathable material with soft hook and loop fasteners for adjustability and comfort; quick clips allow the mask to be removed without changing headgear settings.	Easy-to-adjust soft headgear with EZ Peel tabs.
Interface Description		Small, quiet, and easy to fit and use, the Zest Nasal Mask offers comfort. The design of the Zest Nasal Mask combines the Easy-Clip Silicone Seal and Advanced Air Diffuser with technology that makes Fisher & Paykel Healthcare masks perform well.	One piece of molded silicone. No painful forehead support. No hard plastic frame. Goes below the chin to the top of the nose.	Anatomically shaped face pieces in 5 different sizes that fit into face and under chin. No hard frame, no forehead cushion, patient can wear glasses and read and watch television.	Features a simple design and is approved for multipatient use. Dual-walled silicone cushion is easy to remove and clean. By distributing instead of concentrating pressure, the forehead support minimizes telltale red marks and patient complaints. Mask includes frame with adjustable forehead support, cushion, 4 quick clips, and breathable headgear.	With its patented Auto Seal technology, the EasyLife nasal mask nearly eliminates the need for manual mask adjustments. Its design features a dual-cushion construction—the inner cushion creates an instant, self-adjusting seal, while the outer cushion provides comfortable support.

Company		Philips Respironics	ResMed		Sleepnet Corp	
Product Name		 GoLife for Men	 Mirage FX nasal mask	 Quattro FX full face mask	 IQ Nasal Mask	 Mojo
Website		<a href="http://www.philips.com/GoLife">www.philips.com/GoLife</a>	<a href="http://www.resmed.com">www.resmed.com</a>	<a href="http://www.resmed.com">www.resmed.com</a>	<a href="http://www.sleep-net.com">www.sleep-net.com</a>	<a href="http://www.sleep-net.com">www.sleep-net.com</a>
Price/Warranty		90-day warranty for defects of workmanship and materials	90-day warranty	90-day warranty	30-day comfort guarantee	30-day comfort guarantee
Usage	CPAP	Yes	Yes	Yes	Yes	Yes
	APAP	Yes	Yes	Yes	Yes	Yes
	Bi-level	Yes	Yes	Yes	Yes	Yes
Size Options		P, S, M, L	Standard size fits most. Also available in Mirage FX Wide and Mirage FX for Her.	S, M, L	One size: mold for custom fit	S, M, L. Mold to fit
Latex Free		Yes	Yes	Yes	Yes	Yes
Intentional Leak (ml/s)		4-30 standard liters per minute depending on PAP pressure.	31 l/min @ 10 cm of pressure	37 l/min @ 10 cm of pressure	N/A	N/A
Additional Features		FitPack has all cushions in bag, each wrapped separately.	Squeeze-tab elbow makes it easy to attach and detach the elbow from the frame. Flexible, "butterfly wing" forehead support allows a dynamic fit range without contacting the skin.	Spring Air cushion and spring frame eliminate need for forehead support, creating an open field of vision. Fewer points of adjustment make fitting fast and straightforward.	Softer than skin gel seal. Hand moldable mask frame. Easy, light, and unobstructive headgear.	Softest gel seal. Hand moldable mask frame. Easy-to-use magnetic headgear clips.
Fitting Support Solutions		Stand-alone sizing gauge, Web site video clips.	Mask fitting guides available for download at <a href="http://www.resmed.com">www.resmed.com</a> . Clinical education services, including mask-fitting workshops. For more information, call (800) 424-0737.	Mask fitting guides are available for download at <a href="http://www.resmed.com">www.resmed.com</a> . Clinical education services, including mask-fitting workshops. For more information, call (800) 424-0737.	Step-by-step instructions with diagrams available at company Web site or (800) 742-3646.	Easy-to-use size gauge and step-by-step instructions with diagrams available at <a href="http://www.softgelmask.com">www.softgelmask.com</a> or (800) 222-3727.
Pressure Range		4-30 cm H <sub>2</sub> O	4-20 cm H <sub>2</sub> O	4-30 cm H <sub>2</sub> O	4-20 cm H <sub>2</sub> O	4-20 cm H <sub>2</sub> O
Cleaning Instructions		Mild soapy water without perfumes or lotions. Multiple-patient use lab reprocessing instructions also available.	Wash mask system daily with warm, soapy water. Rinse well. Air dry out of direct sunlight. Do not use scented soaps.	Wash mask system daily with warm, soapy water. Rinse well. Air dry out of direct sunlight. Do not use scented soaps.	For home use, wash with warm water and mild soap.	For home use, wash with warm water and mild soap.
Headgear Description		Simple, preformed headgear. Peel tabs easily adjust in four places for a secure, custom fit.	New, 3D-contoured form-fitting SoftEdge headgear is kind to the face and springs to life, ready to fit and wear.	New contoured headgear combined with the spring frame provides stability, durability, and coverage for comfortable, uninterrupted sleep.	Adjustable Breath-O-Prene headgear.	Adjustable Breath-O-Prene headgear.
Interface Description		GoLife's soft facial contour arms gently hug the face to maintain stability and a secure seal even as the patient moves. Self-adjusting, angled nasal pillows allow a one-step fit.	The Mirage FX is effortless and effective, for comfort through the night. Includes compact frame and soft headgear. With four parts including headgear, wear and care have never been easier.	Full face masks don't have to block a patient's view. Give your patients a mask with a clear field of vision. With a lightweight frame, the Quattro FX offers the coverage and stability of a full face mask in an unobtrusive design.	Hand moldable nasal mask combined with softer than skin gel. This mask can be worn with glasses.	Hand moldable full face mask combined with softer than skin gel. The magnetic clips make attaching, adjusting, and detaching easy.